

Premises Support – Decontamination Service with Medi9 Fogging System

Contact us on 024 7710 3777 or sales@parkmere.com



Background

Parkmere has partnered with a decontamination service provider who has been using the Medi9 fogging system since 2014 and has developed significant expertise in this technique which is 99.9% effective against:

- Spores
- Germs
- Viruses
- Bacteria
- Fungi

Despite the effectiveness against these germ types, unlike alcohols, acids and chlorine derived solutions, Medi9 is non-corrosive and does not cause respiratory or skin irritation.

Also unlike alternative fumigation methods, Medi9 fogging does not require lengthy premises lockdowns.



Where it can be used

Initially only for use within both private and NHS hospitals, with the growing need for high level disinfection services, we are now able to offer Medi9 fogging to commercial clients.

Rapid and effective disinfection services for:

- Schools, colleges and universities
- Public transport, e.g. buses and trains
- Surgeries and clinics, e.g. doctors and dentists
- Essential services, e.g. petrol stations, supermarkets, call centres
- Emergency service vehicles, including fire and ambulance services
- Commercial offices, factories and distribution centres
- Delivery vehicles, including trucks and vans
- Shopping trolleys
- Personnel fogging to protect key workers on workplace exit

Contact us on 024 7710 3777 or sales@parkmere.com



Methodology

The decontamination is carried out using a range of high powered fogging machines which atomise the Medi9 fluid and project it over a large area in a very fine mist, penetrating difficult to clean and frequently missed areas, providing a level of cleanliness that would be difficult and time consuming to replicate with manual cleaning methods.

A targeted approach is used throughout the premises, and the following lists demonstrates just some of the areas and equipment that can be fogged, with the particles penetrating deep into recesses, fixtures and gaps that are not accessible with manual methods...



Areas

- ✓Offices
- ✓Canteens
- ✓Sports halls
- ✓Office seating
- ✓IT suites
- ✓Kitchen areas
- ✓Toilets and washrooms
- ✓Staffrooms
- ✓Sixth form centres
- ✓Nurseries
- ✓Classrooms
- ✓Minibus interiors
- ✓Etc, etc


Equipment

- ✓Desks
- ✓Office chairs and soft seating
- ✓Dining tables
- ✓Classroom furniture
- ✓Catering trolleys
- ✓Sports apparatus
- ✓Blinds
- ✓switches and light fittings
- ✓Bannisters, including underside
- ✓Coat pegs
- ✓Lockers
- ✓Door handles
- ✓Etc, etc



Is it safe to use?

- ✓ Medi9 is safe to use in all environs, there is no need to remove equipment, furniture or soft furnishings.
- ✓ It is carefully calibrated to ensure that the freestanding items, equipment and the fabric of the building will not be damaged through excessive moisture levels.
- ✓ No need to disable the fire alarm or ventilation system.
- ✓ The mist air dries, and the area is safe to access 30 minutes after decontamination has taken place.
- ✓ Medi9 is Alcohol-free, Triclosan-free and Chlorine-free
- ✓ It is non-toxic, non-hazardous, food safe and fragrance free



Data

Medi9 is effective against bacteria, germs, spores, viruses and fungi, including gram positive and gram negative bacteria, enveloped and non enveloped viruses and spore forming bacteria and has been tested to the following standards,

EN1276 – Bactericidal Testing

EN13704 – Sporicidal Testing

EN1650 – Fungicidal Testing

EN14476 – Viricidal Testing

EN1500 – Hygienic Hand Rub Testing



Testing

The Hygiena ATP testing system is used as a verification method to complete the audit cycle for fogging and decontamination services, providing an accurate before and after reading.

Results are maintained in a database and included on the service paperwork, post-visit.



Contact us on 024 7710 3777 or sales@parkmere.com



Action

Contact Parkmere for a quotation – we'll ask for some simple, easy to obtain info so we can provide a quote by rapid return, and if approved, book the service in for an ad-hoc, or a regular scheme.

Phone: 024 7710 3777



Contact us on 024 7710 3777 or sales@parkmere.com



Contact Details:

Santis House
Curriers Close
Charter Avenue Industrial Estate
Coventry
CV4 8AW

Phone: 024 7710 3777
Email: sales@parkmere.com
Web: www.parkmere.com

Contact us on 024 7710 3777 or sales@parkmere.com



Our decontamination service provider has the following accreditations -

ISO9001: 2015, UKAS Certificate 9158

MHRA Registrations, CA006332 and CA011799

Safecontractor and SSIP registered

Cyber Essentials registered

NHS Supply Chain framework participants